Trouble hearing on the phone? Now hear clearly!

CE30[™] Amplified Neckloop

The CE30™ amplified neckloop by Clarity® amplifies up to 30 decibels and is a handsfree solution for T-coil hearing aid users. Included is a 2.5 mm and 3.5mm connector allowing use with most cellular, corded and cordless phones, as well as MP3 players.

See reverse side for more details.





CE30™ Amplified Neckloop

FEATURES

- A hands-free communication solution for hearing aid users
- Designed for T-coil hearing aids
- Included is a 2.5 and 3.5mm connector allowing use with most cellular, corded and cordless phones, as well as MP3 players
- Up to 30 decibels of gain
- · Built in microphone
- Two (2) AAA batteries (included)
- One (I) year limited warranty



Clarity, a Division of Plantronics, Inc. 4289 Bonny Oaks Drive Chattanooga, TN 37406 1-800-426-3738

www.clarityproducts.com

